



Radshield

QUALITY X-RAY PROTECTIVE WEAR & ACCESSORIES

Product Catalogue

Gold Standard
for X-Ray Protective
Wear & Accessories





Radshield

QUALITY X-RAY PROTECTIVE WEAR & ACCESSORIES

Index

Colours & Coverings 1

X-Ray Aprons

Coat Aprons	3
Dental Aprons	4
Dental Bibs	5
Double-Sided Aprons	6
Drop-Off Aprons	7
Half Aprons	8
Maternity Aprons	9
Single-Sided Aprons	10
Vest & Skirt Aprons	11

Apron Carry Bag	12
Apron Cleaning Solution	12
Apron Toggles	12
Arm & Forearm Sleeves	13
Breast Shields	14
Dark Room Safelight	15
Examination Gowns	16
Face Shield	17
Gloves & Mittens	18
Gonad Shields	19
Grabber System	19
L & R Marker Set	20
L&R Stainless Steel Clipon Marker	20
Lead Letters	20
Lead Glasses - Eyewear	21
Lead Blockers	36

Lumbar Support Belt	37
Mobile Height Adj. Shields	38
Ovarian Shields	39
Paediatric Protection	40
Patient Caliper	42
Patient Transfer Board	42
Patient Slider	43
Patient Wrap	43
Positioning Sponges	44
Positioning Bags	44
Protective Cap	45
Protective Shields/Blankets	45
Scoliosis Stole	46
Scrubs	47
Shin Guards	47

Storage


Hangers/Hooks	48
Mobile Apron Trolleys	49
Wall-Mounted Racks	51
Wall-Mounted Peg Racks	52

Straps & Restraints	53
Table Drapes	54
Thyroid Collars Covers	55
Reusable Thyroid Collar Covers	56
Washable Aprons	57
X-Ray Mattresses	57




Colours & Coverings

STANDARD NYLONS




STANDARD VINYLs



STANDARD PVC's



STANDARD TRIMS




RIPSTOP (*Surcharge Applicable*)



Colours & Coverings (Continued)


PRINTS (*Surcharge Applicable*)




LIMITED EDITION (*Surcharge Applicable*)



LIMITED EDITION 2024 FABRICS (*Surcharge Applicable*)



LIMITED EDITION 2025 FABRICS (*Surcharge Applicable*)



X-Ray Aprons

Coat Aprons

Radshield's Coat apron's overlapping front design provides double protection to radiation sensitive sternal and gonadal areas. Worn like a Coat, together with our Support Belt, it is designed for the wearer's total comfort by reducing pressure off the back and shoulders and redistributing the weight.

Available in:

.25mmPb, .35mmPb, .50mmPb or combinations thereof.

Standard Sizes:

Small – XLarge. Larger and Custom Made sizes available upon request.



C1-front



C1-back



C1-belt

Closure:

C1 = Adjustable Hook-and-Loop
Fastener and Support Belt with Quick
Release Buckle

X-Ray Aprons

Dental Aprons

Radshield's Dental Aprons provide your patients with optimum radiation protection, plus they are easy to put on and take off.

Available in: .25mmPb, .35mmPb, .50mmPb or combinations thereof.



D6-front



D6-back

Radshield's Dental Apron without Collar, provides patients with optimum radiation protection. Covers front of body from shoulders to gonadal area.

Standard Sizes:

Small - 54cm Wide x 63cm Long (From Neckline)
Medium - 58cm Wide x 63cm Long (From Neckline)
Large - 62cm Wide x 63cm Long (From Neckline)

Note: Custom made sizes are available upon request.

Closure:

D6 = Adjustable Hook-and-Loop Fastener



D9-front



D9-back

Radshield's Dental Apron complete with Collar, provides patients with optimum radiation protection. Covers front of body from shoulders to gonadal area.

Standard Sizes:

Small - 54cm Wide x 63cm Long (From Neckline)
Medium - 58cm Wide x 63cm Long (From Neckline)
Large - 62cm Wide x 63cm Long (From Neckline)

Note: Custom made sizes are available upon request.

Closure:

D9 = Adjustable Hook-and-Loop Fastener



D17-front



D17-back

Standard Sizes:

Small - 50cm Wide x 50cm Long
Medium - 65cm Wide x 65cm Long
Large - 80cm Wide x 80cm Long
XLarge - 95cm Wide x 95cm Long

Note: Custom made sizes are available upon request.

Closure:

D17 = No Closure



P1-front



P1-back

Radshield's Panorama Apron is required for patient protection when dental x-rays are taken whilst the patient is standing at the panorama machine.

Standard Sizes:

Small – XLarge. Larger and Custom Made sizes available upon request.

Closure:

P1 = Adjustable Hook-and-Loop Fastener

X-Ray Aprons

Dental Aprons (Continued)



DC5-front



DC5-back

Radshield's Dental Cape, provides patients with optimum radiation protection.

Standard Sizes:

Small - 50cm Wide x 50cm Long
Medium - 65cm Wide x 65cm Long
Large - 80cm Wide x 80cm Long
XLarge - 95cm Wide x 95cm Long

Note: Custom made sizes are available upon request.

Closure:

DC5 = Adjustable Hook-and-Loop Fastener



DC6-front



DC6-back

Radshield's Paediatric Dental Cape provides patients with optimum radiation protection.

Available in:

.25mmPb, .35mmPb or .50mmPb.

Standard Size:

Small: 50cm Wide x 50cm Long

Closure:

DC6 = Adjustable Hook-and-Loop Fastener

Note: Smaller and Custom Made sizes available upon request.

Dental Bibs

Radshield's Reversible and Fully Washable Dental Bibs.



D10-front



D10-back

Standard Sizes:

Small - 33cm Wide x 31cm Long (From Neckline)
Medium - 39cm Wide x 41cm Long (From Neckline)
Large - 45cm Wide x 41cm Long (From Neckline)
XLarge - 53cm Wide x 42cm Long (From Neckline)

Note: Dental Bibs can be purchased Individually or as a Set of 4.

Closure:

D10 = Hook-and-Loop at the Back

X-Ray Aprons

Double-Sided Aprons

Radshield's Double-Sided Aprons have been designed for full front and back protection. They fasten at one shoulder for added ease and comfort and both sides are easily and fully adjustable.

Available in:

.25mmPb, .35mmPb, .50mmPb or combinations thereof.

Standard Sizes:

Small – XLarge. Larger and Custom Made sizes available upon request.



D1-front



D1-back



D2-front



D2-back

Closure:

D1 = Adjustable Hook-and-Loop Fastener on Shoulder and on Both Sides

Closure:

D2 = Adjustable Quick Release Buckle on Shoulder and Hook-and-Loop Fastener on Both Sides



D3-front



D3-back

Closure:

D3 = Adjustable Quick Release Buckle on Shoulder and on Both Sides



D18-front



D18-back

Closure:

D18 = Adjustable Quick Release Buckle on Shoulder and in the Front

X-Ray Aprons

Double-Sided Aprons (Continued)



D19-front



D19-back

Radshield's D19 is an apron especially designed to offer complete radiation protection when seated.

Closure:

D19 = Adjustable Quick Release Buckle on Shoulder and in the Front



D20-front



D20-back

Closure:

D20 = Adjustable Quick Release Buckle on Shoulder and Hook-and-Loop Fastener in the Front

Drop-Off Aprons

Radshield's Single and Double-Sided "Drop-Off "Theatre/Surgical Aprons are designed that once the radiological portion of the procedure has been completed, the aprons are released from under the surgical gown and simply fall smoothly to the floor, without compromising sterile fields.

Available in:

.25mmPb, .35mmPb, .50mmPb or combinations thereof.

Standard Sizes:

Small – XLarge. Larger and Custom Made sizes available upon request.



D13-front



D13-back

Closure:

D5 = Adjustable Quick Release Buckle on both Shoulders and Hook-and-Loop Fastener on the Sides

or

D13 = Adjustable Quick Release Buckle on both Shoulders and on the Sides



D16-front



D16-back

Closure:

D4 = Adjustable Hook-and-Loop Fastener in the Front
or

D16 = Adjustable Quick Release Buckle in the Front

X-Ray Aprons

Half Aprons

Radshield's Half Aprons are designed for versatility in patient protection, from infants to adults.

Available in:

.25mmPb, .35mmPb, .50mmPb or combinations thereof.



H1/H6



H18



H18(2)

Standard Sizes:

Small - 35cm (Width) x 30cm (Length)
Medium - 40cm (Width) x 35cm (Length)
Large - 45cm (Width) x 40cm (Length)

Closure:

H1 = Adjustable Quick Release Buckle

or

H6 = Adjustable Hook-and-Loop Fastener

Note: Half Aprons can be purchased Individually or as a Set of 3.

Radshield's "Demi" Half Aprons with flexible Vinyl covered Steel Spring Belt - no Ties or Hook-and-Loop Fastener to adjust. Optional Extra: Demi Rack (**Code: H22**) - see Storage Section.

Standard Sizes:

Small - 20cm (Width) x 25cm (Length)
Medium - 30cm (Width) x 30cm (Length)
Large - 40cm (Width) x 35cm (Length)
XLarge - 50cm (Width) x 40cm (Length)

Closure:

H18 = Flexible Vinyl Steel Spring Belt

Note: Half Aprons can be purchased Individually or as a Set of 4

X-Ray Aprons

Maternity Aprons

Radshield offers Single, Double and Half Aprons in our Maternity range.

Available in:

25mmPb, .35mmPb, .50mmPb or combinations thereof.

Standard Sizes:

Medium – Large. Larger and Custom Made sizes available upon request.



MAT5-front



MAT5-back



MAT4-front



MAT4-back

Radshield's Single-Sided Maternity apron is quick to put on, adjust and take off.

Closure:

MAT5 = Adjustable Quick Release Buckle in the Front

Radshield's Double-Sided Maternity aprons have been designed for full front and back protection. They fasten at one shoulder, for added ease and comfort and both sides are also easily and fully adjustable.

Closure:

MAT2 = Adjustable Hook-and-Loop Fastener on Both Sides

or

MAT4 = Adjustable Quick Release Buckle in the Front



MAT3

Radshield's Maternity Half Apron can be worn under the full apron for double protection in the abdomen/pelvic cavity.

Standard Sizes:

Small - 61cm (Width) x 46cm (Length)

Medium - 61cm (Width) x 54cm (Length)

Large - 61cm (Width) x 61cm (Length)

Closure:

MAT 3 = Adjustable Quick Release Buckle at the Back

X-Ray Aprons

Single-Sided Aprons

Radshield's Single-Sided aprons have always been a popular choice... that's because they're easy to put on, adjust and take off.

Available in:

.25mmPb, .35mmPb, .50mmPb or combinations thereof.

Standard Sizes:

Small – XLarge. Larger and Custom Made sizes available upon request.



S1-front



S1-back



S2-front



S2-back

Closure:

S1 = Adjustable Hook-and-Loop Fastener in the Front

Closure:

S2 = Adjustable Quick Release Buckle in the Front



S3-front



S3-back



S4-front



S4-back

Closure:

S3 = Adjustable Quick Release Buckle at the Back

Closure:

S4 = Adjustable Hook-and-Loop Fastener at the Back

X-Ray Aprons

Single-Sided Aprons (Continued)



S12-front



S12-back



Radshield's S12 apron is an apron especially designed to offer complete radiation protection when seated.

Closure:

S12 = Adjustable Quick Release Buckle in the Front

Vest & Skirt Aprons

Radshield's Vest & Skirt combination has overlapping front panels, thereby providing double protection in critical areas.

Available in:

.25mmPb, .35mmPb, .50mmPb or combinations thereof.

Standard Sizes:

Small – XLarge. Larger and Custom Made sizes available upon request.

Note: Vests and/or Skirts may be purchased as separate items. To order separately, kindly specify size and the piece you prefer. In addition, Vest and/or Skirt sizes may be interchanged (i.e. small Vest with a large Skirt) at no additional charge. Please specify the size for each piece when ordering.



V1-front



V1-back

Closure:

V1 = Adjustable Hook-and-Loop Fastener and Support Belt with Quick Release Buckle

Apron Carry Bag

Radshield's Apron Carry Bag allows you to easily transport your protective apron from one location to another. Constructed of heavy-duty nylon, simply roll your apron around the inner cardboard tube and easily insert it into the barrel-style bag. Complete with adjustable shoulder strap and two handles.



A3

Apron Cleaning Solution

Radshield's X-Ray Apron Cleaning Solution (100% Bio-Degradable) is ideal to disinfect, remove stains and eliminate bad odours. Simply spray the solution onto the apron, rub it into the stain with a soft cloth, rinse with water and allow to dry. Do not submerge apron in water or use an abrasive brush, as it will damage the apron's outer covering.




ACS


Apron Toggles

Radshield offers 2 (two) styles of Apron Toggles, complete with cord.

Note: Sold in Pairs.



T8



T8.1

Arm Sleeve

Radshield's Arm Sleeve provides protection for the arm against scatter radiation. Available in Left, Right or as a Set.

Available in:

.25mmPb, .35mmPb or .50mmPb.




A1



A1-side

HOW TO TAKE THE MEASUREMENTS FOR YOUR **ARM** / **FOREARM SLEEVES**

In order for us to manufacture your Sleeve, we need your measurements as indicated i.e. your bicep **circumference**, as well as the length from the shoulder seam of your shirt/blouse down to your wrist i.e. where you would like it to end ... or your **forearm circumference** measured from below your elbow down to your wrist i.e. where you would like it to end. When measuring, please ensure that your arm is straight and also, that the end of the Sleeve does not cut and/or rub against the top of your hand.



Forearm Sleeve

Radshield's Forearm Sleeve provides protection for the forearm against scatter radiation. Available in Left, Right or as a Set.

Available in:

.25mmPb, .35mmPb or .50mmPb.



F1



F1-side

Breast Shields

Radshield's Adult Breast Shields put the protection where you need it during spinal imaging. Simply adjust the Shields on the Hook-and-Loop Fastener neck strap and position the Shields exactly where you require them.

Available in:

.50mmPb.

Set Comprises of:

2 x Small – 17cm x 20cm

2 x Medium – 20cm x 25cm

2 x Large – 25cm x 30cm

2 x XLarge – 29cm x 35cm



B8

Dark Room Safelight

Radshield's Dark Room Safelight comes complete with Filter and 15 Watt BC Pigmy Globe.
The Safelight is suitable for both green and blue film.

S15 = Safe Light Complete with Filter & Globe

S16 = Safe Light Filter

S17 = Safe Light Globe

Standard Sizes:

Safelight Unit – 27cm Long x 23cm Wide x 11cm Deep

Filter – 25cm Long x 20cm Deep

Note: Can be mounted on Wall or Ceiling.



S15



S15

Examination Gowns

Standard Sizes:

Small – XLarge. Larger and Custom Made sizes available upon request.



E2-front



E2-back



E3-front



E3-back

Closure:

E2 = Ties at the back

Front Wrap Gowns

Closure:

E3 = Ties on the side



E4-front



E4-back

Back Wrap Gowns

Closure:

E4 = Ties on the side

Face Shields

Whether you wear prescription glasses or want to protect your eyes and nose from radiation and liquid splashes, these shields are an indispensable part of your radiation protection attire. The shields extend downwards from a dual adjustable headpiece. Adjust it to the perfect angle for you by simply rotating the dials on the headband for a secure, comfortable fit. The Face Shields cannot be sterilised at high temperatures or with high pressure sterilisation methods, so clean by wiping with a soft cloth and non-abrasive cleaning materials.

Available in: .10mmPb

H20 Weight: 350g

H21 Weight: 395g

Note: Face Shields are supplied with a Storage Bag



H20



H21



Gloves & Mittens



M7

Radshield's M7 Mittens provide complete protection to the hands and are soft, supple and a pleasure to wear.

Available in: .25mmPb, .35mmPb or .50mmPb

Standard Sizes: One Size
Fits All - Sold in Pairs



M8

Radshield's M8 Mittens allow fingers (not thumb) movement, whilst providing complete protection to the thumb and back of the hands. Soft, supple and a pleasure to wear.

Available in: .25mmPb, .35mmPb or .50mmPb

Standard Sizes: One Size
Fits All - Sold in Pairs



M9

Radshield's M9 Mittens allow fingers and thumb movement, whilst providing complete protection to the back of the hands. Soft, supple and a pleasure to wear.

Available in: .25mmPb, .35mmPb or .50mmPb

Standard Sizes: One Size
Fits All - Sold in Pairs



G5

Radshield's G5 Lead Gloves is mainly used in patient positioning but is also ideal for procedures involving high radiation exposure. The glove is soft, flexible, and can easily be wiped clean.

Available in: .50mmPb

Standard Sizes: One Size
Fits All - Sold in Pairs



M16

Radshield's M16 Mittens provide complete protection to the back of the hands, whilst allowing fingers and thumb movement.

Available in: .25mmPb, .35mmPb or .50mmPb

Standard Sizes: One Size
Fits All - Sold in Pairs



M17

Radshield's M17 Mittens provides finger slits for complete finger movement, as well as complete protection to the back of the hands. Adjustable Hook-and-Loop Fastener Strap provides a secure fit.

Available in: .25mmPb, .35mmPb or .50mmPb

Standard Sizes: One Size
Fits All - Sold in Pairs

Gonad Shields



G1

Radshield's Triangular Gonad Shields are specifically for the protection of the genital area.

Available in:
.50mmPb or 1.0mmPb.

Set Comprises of:
1 x Small - 14cm x 14cm x 14cm
1 x Medium - 19cm x 19cm x 19cm
1 x Large - 24cm x 24cm x 24cm

Note: Gonad Shields can be purchased Individually or as a Set of 3.

Optional Extra's:
• Adjustable Hook-and-Loop Fastener Belt (**Code B3**)
• Wall-Mounted Hook-and-Loop Fastener Bar for out of the way storage (**Code VB**)



G18

Radshield's Triangular Gonad Shields Set, complete with Adjustable Webbing Belt are specifically for the protection of the genital area.

Available in:
.50mmPb or 1.0mmPb.

Set Comprises of:
1 x Small - 14cm x 14cm x 14cm
1 x Medium - 19cm x 19cm x 19cm
1 x Large - 24cm x 24cm x 24cm

Note: Gonad Shields can be purchased Individually or as a Set of 3.

Optional Extra:
• Wall-Mounted Hook-and-Loop Fastener Bar for out of the way storage (**Code VB**)



G7

Radshield's Triangular Gonad Shields complete with Ties are specifically for the protection of the genital area.

Available in:
.50mmPb or 1.0mmPb.

Set Comprises of:
1 x Small - 14cm x 14cm x 14cm
1 x Medium - 19cm x 19cm x 19cm
1 x Large - 24cm x 24cm x 24cm

Note: Gonad Shields can be purchased Individually or as a Set of 3.

Optional Extra:
• Wall-Mounted Hook-and-Loop Fastener Bar for out of the way storage (**Code VB**)



VB

Radshield's Hook-and-Loop Fastener Bar (for Wall-Mounting) for the storage of Gonad and Ovarian Shields.



B3

Radshield's fully Adjustable Webbing Belt is used with our Gonad (**Code G1**) and Ovarian (**Code O2**) Shields.

Closure:
B3 = Adjustable Quick Release Buckle

Grabber System

Radshield's Grabber System comes complete with a pair of stainless steel brackets and a wide radiolucent strap to keep patients secure and stable during procedures. Simple to use and easily attaches/detaches to any examination table.



G77

L & R Marker Sets

Radshield's L & R Marker Sets are supplied with a 1 meter cord.

Standard Sizes:

Square

LR3: Small – 15mm x 15mm (2mm thick)

LR4: Large – 25mm x 25mm (2mm thick)

Circular

LR5: Small – 15mm in diameter (2mm thick)

LR7: Medium – 20mm in diameter (2mm thick)

LR6: Large – 25mm in diameter (2mm thick)



L&RMS



L&RMS

L&R Stainless Steel Clip-on Marker Set (Set of 2)

L&R Marker Clip Set



LR8

Lead Letters

Standard Sizes:

LL6: 6mm

Available: A – Z and 0 – 9

LL12: 12mm

Available: A – Z and 0 – 9



LL

Lead Glasses - Eyewear



G19

The **53 Wrap** is constructed of durable nylon that exhibits superior flexibility while minimizing potential breakage. Its full saddle bridge provides maximum comfort and evenly distributes weight over the nose.

A/R Coating

Colour Options: Black, Tortoise

Features

- Spring Hinge temples for enhanced flexibility
- Built-in glass side shields for better field of view and protection
- Optional adjustable nose pads



53 Wrap Eyewear

Characteristics

- Frame Type: Classic
- Front Lens: 0.75mm Pb Glass
- Weight: 84g
- Frame Size: 58-16-135
- Frame Material: Nylon
- Side Shields: 0.50 or 0.75mm Pb Glass
- CE Certified: Yes
- Colour Options: Black, Tortoise
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 62



G31

The **Astroflex Lead goggles** feature molded side shielding and are fitted to the brow for increased splash protection. The goggle's temple length and frame angle are adjustable for a uniform custom fit. The model 70s comes with leaded side shields for necessary lateral x-ray coverage.

A/R Coating

Colour Option: Clear/Blue

Features

- Adjustable temple arms
- Integrated lateral side shields



70s Astroflex Lead Goggles

Characteristics

- Frame Type: Goggles
- Front Lens: 0.75mm Pb Glass
- Weight: 88g
- Frame Size: 59-19-Adjustable
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Strips
- CE Certified: Yes
- Colour Options: Clear/Blue
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 62



G38

The **99 Ultralites** provide a wide area of protection against x-ray. Their high curvature provide adequate lateral protection which eliminates the need for side shields. The 99 comes in 2 sizes and an optional soft rubber nose bridge for additional support and comfort (99R). (Blue only available in regular size and Fuchsia only available in small).

A/R Coating

Colour Options: Brown (Gold), Black, Blue, Silver, Fushia

Sizes: Small, Regular

Features

- Rubberized temples for comfort and grip
- High curvature for maximum field of view and x-ray protection
- Optional rubberized nose bridge (99R)
- Available in 2 size



99 Ultralite Eyewear

Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 63g (regular) | 59.5g (small)
- Frame Size: 61-18-140 (regular) | 56-17-140 (small)
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Strips (optional)
- CE Certified: Yes
- Colour Options: Brown (Gold), Black, Blue, Silver, Fushia
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 64 (regular) 58 (small)



G46

The **Circuit** lead glasses frame features wrap-around lenses which maximize lateral protection without the use of side shields. Rubberized temples and nose pads provide excellent comfort and grip. The Circuits modern style combined with its optimized functionality ensure long lasting quality and protection.

A/R Coating

Colour Options: Silver and Black, Red and Grey, Black and Grey

Features

- Wide Field of view
- Rubberized temples and nose pads for grip and comfort
- Versatile fit for all face types



Circuit Lead Glasses

Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 79g
- Frame Size: 60-18-135
- Frame Material: Polycarbonate
- Side Shields: 0.50mm Pb Strips (optional)
- CE Certified: Yes
- Colour Options: Silver and Black, Red and Grey, Black and Grey
- Rx Options: None
- Min PD Value: Not Applicable

Lead Glasses - Eyewear (Continued)



G48

The **Nike Brazen** lead glasses feature a ventilated nose bridge for grip, comfort, and reduced fogging. The Brazen's frame is constructed of lightweight nylon and offers superior coverage and protection.

A/R Coating

Colour Options: Matte Black/ Blue, Shiny Black

Features

- Ventilated nosebridge
- Lightweight frame
- Superior coverage for protection



Nike Brazen Lead Glasses Characteristics

- Frame Type: (Wraparound)
- Premium Brands – Nike
- Front Lens: 0.75mm Pb Glass
- Weight: 75g
- Frame Size: 60-16-140
- Frame Material: Nylon
- Side Shields: Optional 0.50mm Pb Strips
- CE Certified: Yes
- Color Options: Matte Black/ Blue, Shiny Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 64



G71

The **Wiley-X P-17** lead glasses are low-profile and lightweight, making them a fan favorite. These protective frames boast highly curved lenses for maximum x-ray protection. The P-17 eyewear are ideal for medium to large facial structures.

A/R Coating

Colour Options: Matte Black, Glossy Black, Brown

Features

- Ideal for medium to large face types
- Very durable frame
- Highly curved lenses for maximum x-ray protection



Wiley-X P-17 Lead Glasses Characteristics

- Frame Type: (Wraparound)
- Premium Brands – Wiley-X
- Front Lens: 0.75mm Pb Glass
- Weight: 68g
- Frame Size: 60-18-125
- Frame Material: Polycarbonate
- Side Shields: None
- CE Certified: Yes
- Colour Options: Matte Black, Glossy Black, Brown
- Rx Options: None
- Min Prg PD: Not Applicable



G78

The **Fitguard Alpha** fitover is a versatile, high quality nylon frame. Featuring Flexfit technology, the Fitguard's temples can be easily adjusted for a secure, comfortable fit. The Alpha model is ideal for fitting over wide prescription eyewear while still providing a secure and comfortable fit. With the most balanced weight distribution, the Fitguard Alpha is our most popular fitover model. Blacked out side shields provide a streamlined look as well as 0.50mm lead equivalent protection.

A/R Coating

Colour Options: Matte Black

Features

- 0.50mm lead side shields included
- Ideal for fitting over wide prescription eyewear
- Adjustable temples
- Highly durable frame



Fitguard Alpha Lead Glasses Characteristics

- Frame Type: Fitover
- Front Lens: 0.75mm Pb Glass
- Weight: 91g
- Frame Size: 69-17-145
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Lead Strips
- CE Certified: Yes
- Colour Option: Matte Black
- Rx Options: None
- Min Prg PD: Not Applicable



G79

The **Malibu** lead glasses frame modernizes the classic wayfarer design and is offered in several tasteful color selections. While not one of the most protective styles for lateral and wide-angle x-ray protection, the Malibu is a stylish option for those requiring basic radiation protection at an affordable price.

A/R Coating

Colour Options: Matte Black, Matte Blue, Matte Brown, Shiny Black

Features

- Versatile fit for all face types



Malibu Lead Glasses Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 68.1g
- Frame Size: 54-20-155
- Frame Material: Polycarbonate
- Side Shields: Not available
- CE Certified: Yes
- Colour Options: Matte Blue, Matte Black, Matte Brown, Shiny Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 58

Lead Glasses - Eyewear (Continued)



G81

The **Nike Rabid** lead glasses feature a high-curvature sports frame providing superior wraparound coverage and protection. These protective eyewear boast a lightweight nylon and wide lenses for clarity at all angles during any procedure. Nike's Secure Wrap Technology features grip zones throughout the frame for a better overall fit, comfort, and stability.

A/R Coating

Colour Options: Matte Black/White

Features

- High curvature for wraparound coverage
- Lightweight frame
- Superior grip for better stability and comfort



Nike Rabid Lead Glasses Characteristics

- Frame Type: (Wraparound) Premium Brands – Nike
- Front Lens: 0.75mm Pb Glass
- Weight: 63g
- Frame Size: 63-14-145
- Frame Material: Nylon
- Side Shields: None
- CE Certified: No
- Colour Options: Matte Black/White
- Rx Options: None
- Min PD Value: Not Applicable



G82

A customer favorite, the sporty **Jupiter Squared** is attitude with refinement. The handles of switchblade knives inspired the earstems and their metal bands and rivets. The Matte Black frames of the Flag Collection are highlighted with American flag emblems.

A/R Coating

Colour Option: Matte Black – USA Flag

Features

- Rubberized ear socks for added grip
- and comfort
- American flag on Oakley "O" Logo



Oakley Jupiter Square Lead Glasses Characteristics

- Frame Type: Brand – Oakley
- Front Lens: 0.75mm Pb Glass
- Weight: 74.1g
- Frame Size: 56-18-135
- Frame Material: Durable O Matter™
- Side Shields: 0.50mm Side Strips (optional)
- CE Certified: No
- Colour Option: Matte Black – USA Flag
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 66



G83

The **Razer** Lead Glasses boast a sleek design while being a lightweight and durable frame. These protective eyewear feature a ventilated nose bridge, rubberized temples, and wide-angle lenses which maximize comfort and protection from x-ray. The modern style and use of light-weight materials make the Razer an excellent choice for long-procedures.

A/R Coating

Colour Options: Matte Green, Matte Black, Matte Blue, Matte Grey, Matte Red

Features

- Removable 0.50mm LE side shield insert included
- Rubberized, ventilated nose bridge for maximum comfort
- Rubberized temples for grip
- Lightweight frame



Razer Lead Glasses Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 62.5g
- Frame Size: 58-17-135
- Frame Material: Polycarbonate
- Side Shields: 0.50mm Pb Detachable Side Shields Included
- CE Certified: Yes
- Colour Options: Matte Green, Matte Black, Matte Grey, Matte Blue, Matte Red
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 59



G84

The **Wiley-X Brick** lead glasses offer a flexible rubber temple for support during intense activities whilst the removable Facial Cavity Seal provides maximum splash protection. These protective eyewear boast highly curved lenses for excellent x-ray protection.

A/R Coating

Colour Option: Matte Black

Features

- Flexible, rubberized temple arms
- Highly durable frame
- Removable facial foam cavity seal for added comfort
- Highly curved lenses for maximum x-ray protection



Wiley-X Brick Lead Glasses Characteristics

- Frame Type: (Wraparound) Premium Brands – Wiley-X
- Front Lens: 0.75mm Pb Glass
- Weight: 74g
- Frame Size: 63-18-120
- Frame Material: Polycarbonate
- Side Shields: None
- CE Certified: Yes
- Colour Option: Matte Black
- Rx Options: None
- Min Prg PD: Not Applicable

Lead Glasses - Eyewear (Continued)



G85

The **Centrist** lead glasses are a rectangular-style, safety frame made from high quality nylon and has permanent extended side shields for added x-ray protection. The sleek black with clear insert frame has an adjustable nosepiece that provides comfort and a secure fit.

A/R Coating

Colour Option: Black Crystal

Features

- Extended side shields for additional protection
- Rubberized and adjustable nose pads for grip and comfort



Centrist Lead Glasses Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 88g
- Frame Size: 55-17-140
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Strip plus optional extended side shield
- CE Certified: Yes
- Colour Option: Black Crystal
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 62



G86

The **Wiley-X Boss** lead glasses rest comfortably close to your face while the removable Facial Cavity Seal prevents all small dust and dirt particles from getting into your eyes. These protective eyewear boast highly curved lenses for excellent x-ray protection.

A/R Coating

Colour Options: Kryptek Highlander, Matte Black

Features

- Maximum splash protection
- Highly durable frame
- Removable facial foam cavity seal for added comfort
- Highly curved lenses for maximum x-ray protection



Wiley-X Boss Lead Glasses Characteristics

- Frame Type: (Wraparound)
- Premium Brands – Wiley-X
- Front Lens: 0.75mm Pb Glass
- Weight: 84.5g
- Frame Size: 68-18-125
- Frame Material: Polycarbonate
- Side Shields: None
- CE Certified: Yes
- Colour Options: Kryptek Highlander, Matte Black
- Rx Options: None
- Min Prg PD: Not Applicable



G89

Ray Ban's 4105 wayfarer lead glasses have been constructed into a compact, folding style to be easily portable while still being highly durable.

A/R Coating

Color Option: Matte Black

Features

- Compact folding style for ease of portability
- Durable nylon frame



Ray Ban 4105 Eyewear Characteristics

- Frame Type: (Classic)
- Premium Brands – Ray Ban
- Front Lens: 0.75mm Pb Glass
- Weight: 62g
- Frame Size: 50-23-140
- Frame Material: Nylon
- Side Shields: None
- CE Certified: Yes
- Colour Option: Matte Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 65



G90

The **9941 Ultralite** lead glasses come standard with a rubberized nose bridge. These protective eyewear are designed to offer comfort and provide protection for the entire eye area without the need of additional lead side shields. The 9941 glasses are a popular frame and are ideal for those who have a high power prescription. Available in 2 sizes and optional adjustable nose pads for regular sized frame (9941A).

A/R Coating

Colour Options: Red, Silver, Black, Blue

Features

- Rubberized temples for comfort and grip
- Ideal for high prescriptions
- Available in regular and small sizes
- Adjustable nose pads available



9941 Ultralite Lead Glasses

Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 64g (regular) | 60g (small)
- Frame Size: 55-22-130 (regular) | 53-22-130 (small)
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Strips optional
- CE Certified: Yes
- Color Options: Red, Silver, Black, Blue
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 64 (regular) | 58 (small)

Lead Glasses - Eyewear (Continued)



G91

The **98 Superlite** is an ultra light-weight frame designed for comfort. Its high curvature provides enhanced lateral/ side protection against x-ray. Rubberized temples provide a secure, comfortable fit.

A/R Coating

Colour Option: Clear Brown

Features

- Rubberized temples and nose pads for comfort and grip
- High curvature for maximum field of view and x-ray protection
- Lightweight frame



98 Superlite Eyewear

Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 63g
- Frame Size: 63-15-120
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Strips (optional)
- CE Certified: Yes
- Colour Option: Clear Brown
- Rx Options: None
- Min PD Value: Not Applicable



G92

The **Airborne** lead glasses come fully loaded with an Asian fit nose bridge, removable facial cavity, and curved lenses for optimal lateral x-ray protection.

A/R Coating

Colour Option: Clear/Black

Features

- Removable facial cavity foam insert for added comfort
- High curvature for maximum x-ray protection and wide view of view
- Rubberized Asian fit nose bridge for grip and comfort



Airborne Lead Glasses

Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 74.7g
- Frame Size: 57-21-145
- Frame Material: Polycarbonate
- Side Shields: 0.50mm Pb Strips (optional)
- CE Certified: Yes
- Colour Option: Clear/Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 62



G93

The **Mako** lead glasses frame is boat captain approved with its unique patterns and high curve lenses which maximize radiation protection coverage. Integrated side shields equipped with lead strip side shields provide essential lateral x-ray protection.

A/R Coating

Colour Option: Grey/Green

Features

- High curvature for maximum field of view and x-ray protection
- Rubberized temples and nose pads for grip and comfort
- Integrated 0.50mm Pb lateral side shields



Mako Lead Glasses

Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 69.3g
- Frame Size: 65-15-130
- Frame Material: Polycarbonate
- Side Shields: 0.50mm Pb Strips included
- CE Certified: Yes
- Colour Option: Grey/Green
- Rx Options: None
- Min PD Value: Not Applicable



G94

Deep stylistic curves and lenses combine style with lateral x-ray radiation protection. The **Onyx** frame features a ventilated nose bridge providing grip, reduced fogging, and comfort.

A/R Coating

Colour Options: Matte Black, Shiny White

Features

- High curvature for maximum field of view and x-ray protection
- Rubberized, ventilated nose bridge for maximum comfort and grip



Onyx Lead Glasses

Characteristics

- Frame Type: Wraparound
- Front Lens: 0.75mm Pb Glass
- Weight: 71.8g
- Frame Size: 68-19-135
- Frame Material: Polycarbonate
- Side Shields: None
- CE Certified: Yes
- Color Options: Matte Black, Shiny White
- Rx Options: None
- Min PD Value: Not Applicable

Lead Glasses - Eyewear (Continued)



G95

The **Modern Metal Lead Glasses** feature a light-weight, thin wire frame with contoured 0.50mm LE side shields for maximum coverage. Non-slip, silicone nose pads offer easy adjustability and allow the Modern Metal frame to fit a wide variety of face types.

A/R Coating

Colour Options: Gold, Light Gunmetal

Features

- Included 0.50mm LE side shields
- Adjustable silicone nose pads
- Versatile fit for all face types



Modern Metal Lead Glasses Characteristics

- Frame Type: Classic Metals
- Front Lens: 0.75mm Pb Glass
- Weight: 73g
- Frame Size: 48-20-135
- Frame Material: Optical Steel
- Side Shields: 0.50mm LE side shields included
- CE Certified: Yes
- Colour Options: Gold, Light Gunmetal
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 58



G96

The **Metal Aviator Lead Glasses** are a stylish frame that come standard with contoured 0.50mm LE side shields for maximum lateral coverage. Non-slip, silicone nose pads offer easy adjustability and allow the Modern Metal frame to fit a wide variety of face types.

A/R Coating

Colour Option: Gold

Features

- Included 0.50mm LE side shields
- Adjustable silicone nose pads
- Versatile fit for all face types



Metal Aviator Lead Glasses Characteristics

- Frame Type: Classic Metals
- Front Lens: 0.75mm Pb Glass
- Weight: 82g
- Frame Size: 54-17-150
- Frame Material: Optical Steel
- Side Shields: 0.50mm LE side shields included
- CE Certified: Yes
- Colour Option: Gold
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 64



G97

The **Petite Metal lead glasses** are made a smaller, lightweight frame made of optical steel with bronze electroplate. This durable frame comes standard with contoured 0.50mm LE side shields for maximum lateral coverage. Non-slip, silicone nose pads offer easy adjustability and allow the Modern Metal frame to fit a wide variety of face types.

A/R Coating

Colour Options: Gold, Light Gunmetal

Features

- Included 0.50mm LE side shields
- Adjustable silicone nose pads
- Petite, lightweight frame



Petite Metal Lead Glasses Characteristics

- Frame Type: Classic Metals
- Front Lens: 0.75mm Pb Glass
- Weight: 63g
- Frame Size: 52-22-150
- Frame Material: Optical Steel
- Side Shields: 0.50mm LE side shields included
- CE Certified: Yes
- Color Options: Gold, Light Gunmetal
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 58



G98

The **Classic Metal lead glasses** are a lightweight frame made of optical steel with bronze electroplate. This durable frame comes standard with contoured 0.50mm LE side shields for maximum lateral coverage. Non-slip, silicone nose pads offer easy adjustability and allow the Modern Metal frame to fit a wide variety of face types.

A/R Coating

Colour options: Nano Blue, Gold

Features

- Included 0.50mm LE side shields
- Adjustable silicone nose pads
- Lightweight frame



Classic Metal Lead Glasses Characteristics

- Frame Type: Classic Metals
- Front Lens: 0.75mm Pb Glass
- Weight: 64g
- Frame Size: 54-19-140
- Frame Material: Optical Steel
- Side Shields: 0.50mm LE side shields included
- CE Certified: Yes
- Colour Options: Nano Blue, Gold
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 58

Lead Glasses - Eyewear (Continued)



G99

The model **553** is a metal frame featuring custom fitted 0.50mm lead strip side shields that create a seal to the front lenses for complete lateral coverage. Non-slip, silicone nose pads provide easy adjustability and allow the frame to fit a wide variety of face types.

A/R Coating

Colour Option: Silver

Features

- Included 0.50mm LE side shields
- Adjustable silicone nose pads
- Versatile fit for all face types



553s Metalite Eyewear Characteristics

- Frame Type: Classic Metals
- Front Lens: 0.75mm Pb Glass
- Weight: 62g
- Frame Size: 50-18-140
- Frame Material: Optical Steel
- Side Shields: 0.50mm Pb Strips Included
- CE Certified: Yes
- Colour Option: Silver
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 58



G100

The model **557** is a metal frame featuring custom fitted 0.50mm lead strip side shields that create a seal to the front lenses for complete lateral coverage. Non-slip, silicone nose pads provide easy adjustability and allow the frame to fit a wide variety of face types.

A/R Coating

Colour Option: Gold

Features

- Included 0.50mm LE side shields
- Adjustable silicone nose pads
- Versatile fit for all face types



557S Metalite Eyewear Characteristics

- Frame Type: Classic Metals
- Front Lens: 0.75mm Pb Glass
- Weight: 63g
- Frame Size: 48-20-135
- Frame Material: Optical Steel
- Side Shields: 0.50mm Pb Strips Included
- CE Certified: Yes
- Colour Option: Gold
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 58



G101

The model **661/662S** is a metal frame featuring custom fitted 0.50mm lead strip side shields that create a seal to the front lenses for complete lateral coverage. Non-slip, silicone nose pads provide easy adjustability and allow the frame to fit a wide variety of face types.

A/R Coating

Colour Option: Gold

Features

- Included 0.50mm LE side shields
- Adjustable silicone nose pads
- Versatile fit for all face types



661 Metalite Eyewear Characteristics

- Frame Type: Classic Metals
- Front Lens: 0.75mm Pb Glass
- Weight: 54g (no side shields) | 70g (with side shield)
- Frame Size: 52-18-144
- Frame Material: Optical Steel
- Side Shields: 0.50mm Pb Strips Optional (Model 662S)
- CE Certified: Yes
- Colour Option: Gold
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 62



G103


Designed to fit over prescription eyewear, the **90 fitover** lead glasses are an excellent choice for those requiring protection against x-ray and liquid splashes. The large, clear side shields protect the user from lateral exposure while eliminating the "tunnel" effects of look-alikes.

A/R Coating

Colour Option: Matte Black

Features

- 0.50mm glass lead side shields included
- Wide angle field of view
- Large glass side shields help reduce "tunnel effect"



90 Fitover Eyewear Characteristics

- Frame Type: Fitover
- Front Lens: 0.75mm Pb Glass
- Weight: 112g
- Frame Size: 66-14-120
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Lead Glass or 0.75mm Pb Glass
- CE Certified: Yes
- Colour Options: Matte Black
- Rx Options: None
- Min Prg PD: Not Applicable

Lead Glasses - Eyewear (Continued)



G105

The **Fitguard Beta** fitover is a versatile, high quality nylon frame. Featuring Flexfit technology, the Fitguard lead glasses' temples can be benched (adjusted) easily for a secure, comfortable fit. The Beta model is ideal for fitting over standard prescription eyewear with square or rounded lenses. Blacked out side shields provide a streamlined look as well as 0.50mm lead equivalent protection.

A/R Coating

Colour Option: Matte Black

Features

- 0.50mm lead side shields Included
- Ideal for fitting over eyewear with square or rounded lenses
- Adjustable temples
- Highly durable frame



Fitguard Beta Lead Glasses Characteristics

- Frame Type: Fitover
- Front Lens: 0.75mm Pb Glass
- Weight: 108g
- Frame Size: 64-18-150
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Lead Strips
- CE Certified: Yes
- Colour Option: Matte Black
- Rx Options: None
- Min Prg PD: Not Applicable



G106

Similar to the Proguard LeadR's frame, the **Comet** Lead goggles are an affordable option providing excellent coverage against x-ray. The Comet's unique side shield configuration allows for lead shielding to be placed on the frame exterior and wrap around the entire lens area, providing a complete seal from scatter radiation.

A/R Coating

Colour Options: Blue, Pink, Grey/Orange

Features

- Rubberized nose bridge and temples for comfort and fit
- Extended lateral side shields
- Foam top bar for comfort
- Wide field of view



Comet Eyewear Characteristics

- Frame Type: Goggles
- Front Lens: 0.75mm Pb Glass
- Weight: 78.3g
- Frame Size: 59-19-114
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Strips, Optional extended side shields
- CE Certified: Yes
- Colour Options: Blue, Pink, Grey/Orange
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 64



G107

The LeadR Lead Goggles feature a soft, ventilated nose bridge with rubberized temples for premium comfort. Each **LeadR** eyewear frame comes with the added protection of 0.50mm lead side shields and the ability to extend that side protection into an additional panel.

A/R Coating

Colour Options: Black, Grey/Pink, Grey/Yellow, White/Purple

Features

- Extended lateral side shields
- Rubberized temples for comfort and fit
- Rubberized, ventilated nose bridge
- Foam top bar for added comfort



Proguard LeadR Lead Goggles Characteristics

- Frame Type: Goggles
- Front Lens: 0.75mm Pb Glass
- Weight: 112g
- Frame Size: 66-14-120
- Frame Material: Nylon
- Side Shields: 0.50mm Pb Strips
- CE Certified: Yes
- Colour Options: Black, Gray/Pink, Green/Yellow, White/Purple
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 62



G108

The **Oakley Five Squared** lead glasses use Condensed Cranial Geometry for a tapered architecture that feels as good as it looks. Dimensional reliefs, metal icons and sculpturally integrated hinges take the whole idea of street style and amp it with attitude.

A/R Coating

Colour Options: Polished Black, Grey Smoke

Features

- Lightweight O Matter™ frame



Oakley Five Squared Lead Glasses Characteristics

- Frame Type: (Wraparound) Premium Brands – Oakley
- Front Lens: 0.75mm Pb Glass
- Weight: 71g
- Frame Size: 54-20-135
- Frame Material: O Matter™
- Side Shields: None
- CE Certified: Yes
- Colour Options: Polished Black, Grey Smoke
- Rx Options: None
- Min Prg PD: Not Applicable

Lead Glasses - Eyewear (Continued)



G109

The **Abaco Dockside** lead glasses have a bold square style in an impact-resistant, polycarbonate frame featuring premium metal hinges. These fashion-forward and protective frames are both durable and lightweight.

A/R Coating

Colour Options: Black, Crystal Clear

Features

- Unisex frame fits all face types
- Lightweight frame



Abaco Dockside Eyewear Characteristics

- Frame Type: (Classic) Premium Brands – Abaco
- Front Lens: 0.75mm Pb Glass
- Weight: 72g
- Frame Size: 58-19-135
- Frame Material: Polycarbonate
- Side Shields: None
- CE Certified: Yes
- Colours Options: Black, Crystal Clear
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 60



G110

The Abaco Jesse lead glasses are unisex, universally flattering frames that are made of lightweight polycarbonate. These protective eyewear feature characteristic rivets at the temples and extremely flexible spring hinges for a snug, comfortable fit.

A/R Coating

Colour Option: Matte Black

Features

- Flexible spring hinge temples
- Lightweight frame
- Unisex frame fits all face types



Abaco Jesse Eyewear Characteristics

- Frame Type: (Classic) Premium Brands – Abaco
- Front Lens: 0.75mm Pb Glass
- Weight: 65g
- Frame Size: 56-18-135
- Frame Material: Polycarbonate
- Side Shields: 0.50mm Pb Strips
- CE Certified: Yes
- Colour Option: Matte Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 60



G111

Made of impact-resistant, lightweight polycarbonate, the Abaco Kai Lead glasses feature premium metal hinges and are molded around the ears to provide a comfortable fit. These unisex, protective frames complement and fit virtually all faces from small to large.

A/R Coating

Colour Options: Matte Black, Crystal Clear

Features

- Unisex frame fits all face types
- Lightweight frame



Abaco Kai Eyewear Characteristics

- Frame Type: (Classic) Premium Brands – Abaco
- Front Lens: 0.75mm Pb Glass
- Weight: 65g
- Frame Size: 56-18-135
- Frame Material: Polycarbonate
- Side Shields: 0.50mm Pb Strips
- CE Certified: Yes
- Colour Options: Matte Black, Crystal Clear
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 60



G112

The **Abaco Tiki** is a classic-style frame made of a wood-textured, durable and lightweight polycarbonate. These lead glasses feature characteristic rivets at the temples, standard hinges, and are molded around the ears for added comfort. The Tiki protective eyewear are unisex and universally flattering to fit all faces from small to large.

A/R Coating

Colour Option: Dark Wood

Features

- Unisex frame fits all face types
- Lightweight frame



Abaco Tiki Lead Glasses Characteristics

- Frame Type: (Classic) Premium Brands – Abaco
- Front Lens: 0.75mm Pb Glass
- Weight: 70g
- Frame Size: 54-21-140
- Frame Material: Polycarbonate
- Side Shields: 0.50mm Pb Strips, Optional Extended Side Shields
- CE Certified: Yes
- Colour Options: Dark Wood
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 60

Lead Glasses - Eyewear (Continued)



G113

The classic Abaco Waikiki lead glasses feature a lightweight polycarbonate frame to provide durability and a comfortable fit. The Abaco leaded eyewear feature characteristic rivet details on the front, premium metal hinges, and they are molded around the ears for added comfort. Waikiki glasses look good on virtually everyone.

A/R Coating

Colour Options: Tortoise, Gloss Black

Features

Unisex frame fits all face types
Lightweight frame



Abaco Waikiki Lead Glasses Characteristics

- Frame Type: (Classic) Premium Brands – Abaco
- Front Lens: 0.75mm Pb Glass
- Weight: 68g
- Frame Size: 54-21-140
- Frame Material: Polycarbonate
- Side Shields: N/A
- CE Certified: Yes
- Colour Options: Tortoise, Gloss Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 60



G114

The Nike Brazen Fury features rubberized nose pads and temple tips for added grip and comfort. With super-comfortable features in a lightweight frame, they are an ideal pair for active and everyday wear.

A/R Coating

Colour Options: Matte Black with Green, Matte Gray with Blue

Features

- Rubberized nose pads and temple tips
- Lightweight and durable plant-based frame (bioplastics)



Nike Brazen Fury Glasses Characteristics

- Frame Type: Brand – Nike
- Front Lens: 0.75mm Pb Glass
- Weight: –
- Frame Size: 60-17-130
- Frame Material: Bioplastics
- Side Shields: 0.50mm Pb Strips (optional)
- CE Certified: No
- Color Options: Matte Black with Green, Matte Gray with Blue
- Rx Options: Single Vision, Bifocal, Progressive
- Min PRG PD: 64



G115

Utilizing Nike's unique nose bridge design, the Maverick lead glasses have reduced fogging and maximized comfort during long fluoroscopic procedures.

A/R Coating

Colour Option: Matte Black

Features

- Unique nose bridge design for fog reduction
- Lightweight frame



Nike Maverick Lead Glasses Characteristics

- Frame Type: Premium Brands – Nike
- Front Lens: 0.75mm Pb Glass
- Weight: 70.3g
- Frame Size: 59-15-145
- Frame Material: Nylon
- Side Shields: None
- CE Certified: No
- Colour Option: Matte Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 64



G116

The Nike Premier lead glasses offer a diverse selection of colors and tasteful design elements. The Premier protective eyewear feature large rectangular lenses to provide adequate x-ray protection and coverage.

A/R Coating

Colour Options: Matte Tortoise, Matte Anthracite, Matte Black

Features

- Large Rectangular Lenses
- Lightweight Frame



Nike Premier Lead Glasses Characteristics

- Frame Type: (Classic) Premium Brands – Nike
- Front Lens: 0.75mm Pb Glass
- Weight: 67.8g
- Frame Size: 60-13-135
- Frame Material: Nylon
- Side Shields: None
- CE Certified: No
- Colour Options: Matte Tortoise, Matte Anthracite, Matte Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min PD Value: 64

Lead Glasses - Eyewear (Continued)



G117

The **Liberty Maxx's** hingeless temples are ultra-sturdy while the head strap keeps these glasses fitting securely during long procedures. Soft padding inside the temples and bridge ensures maximum comfort and the best possible fit.

A/R Coating

Colour Optionst: Shiny Black, Shiny Blue, Shiny Red, Shiny Silver

Features

- Available in 2 sizes
- Hingeless temples for a secure fit
- Can accommodate wide range of prescriptions
- Versatile fit for all face types (unisex)
- Rubberized nose bridge and temples for added comfort and fit



Liberty Maxx30 Lead Glasses Characteristics

- Frame Type: Brand – Liberty
- Front Lens: 0.75mm Pb Glass
- Weight: 60g (small) | 69g (reg.)
- Frame Size: 53-17-130 (small) 55-20-130 (reg.)
- Frame Material: Polycarbonate
- Side Shields: 0.50mm Side Strips (optional)
- CE Certified: No
- Colour Option: Shiny Black, Shiny Blue, Shiny Red, Shiny Silver
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 56



G118

The **Liberty Maxx's** head strap keeps these glasses fitting securely during long procedures. Soft padding inside the temples and bridge ensures maximum comfort and the best possible fit.

A/R Coating

Colour Options: Shiny Blue, Shiny Black, Shiny Silver

Features

- Available in 2 sizes
- Shortened temples with a fitted head strap
- Can accommodate wide range of prescriptions
- Versatile fit for all face types (unisex)
- Rubberized nose bridge for comfort



Liberty Maxx31 Lead Glasses Characteristics

- Frame Type: Brand – Liberty
- Front Lens: 0.75mm Pb Glass
- Weight: 60g (small) | 64g (reg.)
- Frame Size: 53-17-130 (small) 55-20-130 (reg.)
- Frame Material: Polycarbonate
- Side Shields: 0.50mm Side Strips (optional)
- CE Certified: Yes
- Colour Option: Shiny Blue, Shiny Black, Shiny Silver
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 56



G119

The **Oakley Crosslink Zero** lead glasses have fixed temples with a three-Point Fit to ensure lenses are aligned for visual correction. Wire cores hidden inside the temples of this lightweight frame offer adjustability for a perfect fit.

A/R Coating

Colour Options: Matte Black, Polished Grey Shadow

Features

- Rubberized nosepads and ear socks for comfort
- Lightweight O Matter™ frame



Oakley Crosslink Zero Lead Glasses Characteristics

- Frame Type: (Wraparound) Premium Brands – Oakley
- Front Lens: 0.75mm Pb Glass
- Weight: 61.6g
- Frame Size: 54-20-135
- Frame Material: O Matter™
- Side Shields: None
- CE Certified: Yes
- Colour Options: Matte Black, Polished Grey Shadow
- Rx Options: Single Vision, Min PD Value: 56



G120

The **Oakley Five Squared** lead glasses use Condensed Cranial Geometry for a tapered architecture that feels as good as it looks. Dimensional reliefs, metal icons and sculpturally integrated hinges take the whole idea of street style and amp it with attitude.

A/R Coating

Colour Options: Polished Black, Grey Smoke

Features

- Lightweight O Matter™ frame



Oakley Five Squared Lead Glasses Characteristics

- Frame Type: (Wraparound) Premium Brands – Oakley
- Front Lens: 0.75mm Pb Glass
- Weight: 71g
- Frame Size: 54-20-135
- Frame Material: O Matter™
- Side Shields: None
- CE Certified: Yes
- Colour Options: Polished Black, Grey Smoke
- Rx Options: None
- Min Prg PD: Not Applicable

Lead Glasses - Eyewear (Continued)



G121

The **Oakley Holbrook** lead glasses combine traditional design with a modern flare.

A/R Coating

Colour Options: Polished Clear, Satin Black, Matte Brown Tortoise

Features

- Lightweight O Matter™ frame



Oakley Holbrook Lead Glasses Characteristics

- Frame Type: (Classic) Premium Brands – Oakley
- Front Lens: 0.75mm Pb Glass
- Weight: 66.5g
- Frame Size: 56-18-135
- Frame Material: O Matter™
- Side Shields: None
- CE Certified: No
- Colour Options: Polished Clear, Satin Black, Matte Brown Tortoise
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 60



G122

Smooth details and streamlined design provide a timeless look with the **Oakley Litebeam** lead glasses. Composed of lightweight O Matter™ material, the Litebeam protective eyewear feature adjustable temples and come with several nosebridge sizes to ensure maximum comfort.

A/R Coating

Colour Option: Satin Black

Features

- Adjustable temples for a perfect fit
- Lightweight O Matter™ frame



Oakley Litebeam Eyewear Characteristics

- Frame Type: (Classic) Premium Brands – Oakley
- Front Lens: 0.75mm Pb Glass
- Weight: 59.6g
- Frame Size: 55-18-140
- Frame Material: O Matter™
- Side Shields: None
- CE Certified: No
- Colour Option: Satin Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 60



G123

The **Oakley Straightlink** lead glasses feature rubberized earsocks and nose pads for a no-slip grip with perspiration, making this frame great for long procedures. These protective glasses have a Three-Point Fit which holds lenses in precise optical alignment while eliminating pressure points. The Straightlink's O Matter™ stress-resistant frame material is both lightweight and durable for all-day comfort and protection.

A/R Coating

Colour Option: Matte Black

Features

- Rubberized earsocks and nose pads for comfort and grip
- Optimized weight distribution and alignment
- Lightweight O Matter™ frame



Oakley Straightlink Lead Glasses Characteristics

- Frame Type: (Wraparound) Premium Brands – Oakley
- Front Lens: 0.75mm Pb Glass
- Weight: 77.2g
- Frame Size: 61-19-130
- Frame Material: O Matter™
- Side Shields: None
- CE Certified: Yes
- Colour Option: Matte Black
- Rx Options: None
- Min Prg PD: Not Applicable



G124

Designed for small to medium faces, **Oakley's Trillbe X** lead glasses have a rounded front with a sculptured keyhole bridge. These protective eyewear's architecture take vintage style into modern territory.

A/R Coating

Colour Options: Matte Black, Polished Clear

Features

- Designed for small to medium faces
- Lightweight frame



Oakley Trillbe X Lead Glasses Characteristics

- Frame Type: (Classic) Premium Brands – Oakley
- Front Lens: 0.75mm Pb Glass
- Weight: 64g
- Frame Size: 52-18-140
- Frame Material: O Matter™
- Side Shields: None
- CE Certified: Yes
- Colour Options: Matte Black, Polished Clear
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 58

Lead Glasses - Eyewear (Continued)



G125

Ray Ban's Daddy-O lead glasses are a sporty frame that boast optimal coverage for maximum x-ray protection.

A/R Coating

Colour Option: Matte Black

Features

- Excellent coverage for maximum protection
- Lightweight frame



Ray Ban Daddy-O Lead Glasses Characteristics

- Frame Type: (Wraparound) Premium Brands – Ray Ban
- Front Lens: 0.75mm Pb Glass
- Weight: 61g
- Frame Size: 59-17-125
- Frame Material: Nylon
- Side Shields: Optional 0.50mm Pb Strips
- CE Certified: Yes
- Colour Option: Matte Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 65



G126

Ray Ban's 3534 lead glasses boast a sporty, masculine frame perfect for everyday use. Brushed aluminum temples, adjustable nose bridge, and rubberized ear socks ensure your eyewear are perfectly balanced and have a personalized fit.

A/R Coating

Colour Options: Gunmetal, Black

Features

- Adjustable silicone nose pads for personalized fit
- Rubberized ear socks for comfort



Ray Ban 3534 Lead Glasses Characteristics

- Frame Type: Premium Brands – Ray Ban
- Front Lens: 0.75mm Pb Glass
- Weight: 63g
- Frame Size: 59-20-135
- Frame Material: Metal
- Side Shields: None
- CE Certified: Yes
- Colour Options: Gunmetal, Black
- Rx Options: None
- Min Prg PD: Not Applicable



G127

Inspired by the iconic Wayfarer sunglasses, **Ray Ban's 4165** Justin lead glasses are redesigned with slightly larger rectangle lenses and a rubberized frame finish to really give it that extra edge.

A/R Coating

Colour Option: Matte Black

Features

- Rubberized frame finish for no-slip grip and comfort
- Durable Nylon Frame



Ray Ban 4165F Lead Glasses Characteristics

- Frame Type: (Classic) Premium Brands – Ray Ban
- Front Lens: 0.75mm Pb Glass
- Weight: 73g
- Frame Size: 55-18-145
- Frame Material: Nylon
- Side Shields: None
- CE Certified: Yes
- Color Options: Matte Black
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 65



G128

Designed to fit all face types, **Ray Ban's 4181** lead glasses feature a larger plastic frame that is both comfortable and stylish.

A/R Coating

Colour Option: Tortoise

Features

- Versatile fit for all face types (unisex)



Ray Ban 4181 Eyewear Characteristics

- Frame Type: (Classic) Premium Brands – Ray Ban
- Front Lens: 0.75mm Pb Glass
- Weight: 71.4g
- Frame Size: 57-16-145
- Frame Material: Injected Plastic
- Side Shields: None
- CE Certified: No
- Colour Option: Tortoise
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 65

Lead Glasses - Eyewear (Continued)



G129

A sharper more self assured Ray Ban style, the model **7047** is full of confidence and style. Composed of lightweight materials, the 7047 provides comfortable, long lasting protection against x-ray.

A/R Coating

Colour Options: Black, Blue, Brown

Features

- Wide fit temples for a comfortable fit



Ray Ban 7047 Lead Glasses Characteristics

- Frame Type: Brand – Ray Ban
- Front Lens: 0.75mm Pb Glass
- Weight: 63g
- Frame Size: 54-17-140
- Frame Material: Injected Plastic
- Side Shields: N/A
- CE Certified: No
- Colour Options: Black, Blue, Brown
- Rx Options: N/A



G130

Designed to withstand high impact, the **Wiley-X Airrage** lead glasses feature a facial foam cavity seal and high curved lenses for maximum comfort and protection against x-ray.

A/R Coating

Colour Option: Matte Black

Features

- Maximum splash protection
- Highly durable frame
- Facial foam cavity for added comfort
- Highly curved lenses for maximum x-ray protection



Wiley-X Airrage Lead Glasses Characteristics

- Frame Type: (Wraparound) Premium Brands – Wiley-X
- Front Lens: 0.75mm Pb Glass
- Weight: 67.5g
- Frame Size: 61-18-124
- Frame Material: Polycarbonate
- Side Shields: None
- CE Certified: No
- Colour Option: Matte Black
- Min Prg PD: Not Applicable



G131

Lightweight and comfortable, the **Wiley-X Contour** lead glasses feature removable side shield inserts to ensure full lateral protection from x-ray.

A/R Coating

Colour Options: Matte Black, Gloss Black/ Brown Stripe

Features

- Removable 0.50mm side shield inserts
- Durable, lightweight frame



Wiley-X Contour Lead Glasses Characteristics

- Frame Type: (Classic) Premium Brands – Wiley-X
- Front Lens: 0.75mm Pb Glass
- Weight: 62.3g
- Frame Size: 54-17-140
- Frame Material: Nylon
- Side Shields: Removable 0.50mm Pb side shield inserts
- CE Certified: Yes
- Colour Options: Matte Black, Gloss Black/ Brown Stripe
- Rx Options: Single Vision, Bifocal, Progressive
- Min Prg PD: 60



G132

Radshield's **Tugga Baby** pediatric glasses provide essential protection during medical fluoroscopy and feature a universal fit with its adjustable headband. These eyewear also come with a matching, Tugga baby eyewear case and a microfiber cleaning cloth.

Colour Options: Pink, Green, Blue

Features

- Adjustable headband for universal fit
- Lightweight frame for small faces



Tugga Baby Lead Glasses Characteristics

- Frame Type: Goggles
- Front Lens: 0.75mm Pb Glass
- Weight: 42g
- Frame Size: 46-18-Adjustable
- Frame Material: Polycarbonate
- Side Shields: Not available
- CE Certified: Yes
- Colour Options: Pink, Green, Blue
- Rx Options: None

Lead Glasses - Eyewear (Continued)



H20

The **Panoramic Acrylic Face Shield** utilizes a lead impregnated acrylic shield suspended from a dual adjustable headpiece. This lead face shield provides additional x-ray protection due to its larger lead impregnated shield. Because of the multi-adjustable headpiece, you can angle the mask to protect from radiation in many different positions as well as adjust the size to comfortably fit your head.

Mask should be worn with leaded eyewear protection.

Features

- Maximum x-ray protection
- Dual-adjustable head piece for guaranteed fit



Panoramic Acrylic Face Shield Characteristics

- Frame Type: Face Shield
- Lens Material: >0.1mm Pb LE Acrylic
- Weight: 350g
- Frame Size: Adjustable 20" – 25" circumference
- CE Certified: Yes
- Rx Options: None



H21

The **Full Acrylic Face Shield** utilizes a lead impregnated acrylic shield suspended from a dual adjustable headpiece. This lead face shield provides additional x-ray protection due to its larger lead impregnated shield. Because of the multi-adjustable headpiece, you can angle the mask to protect from radiation in many different positions as well as adjust the size to comfortably fit your head.

Mask should be worn with leaded eyewear protection.

Colour Option: Black/Clear

Features


- Large protective field
- Dual-adjustable head piece for guaranteed fit



Full Acrylic Face Shield Characteristics

- Frame Type: Face Shield
- Lens Material: >0.1mm Pb LE Acrylic
- Weight: 450g
- Frame Size: Adjustable 20" – 25" circumference
- CE Certified: Yes
- Colour Option: Black/Clear
- Rx Options: None

Lead Blockers



B9

Radshield's Self-Coated Lead Blockers have a smooth non-absorbent vinyl surface that can be easily cleaned.

Available in:
.50mmPb, 1.0mmPb.

Standard Sizes:
24cm x 10cm
30cm x 12cm
43cm x 18cm

Note:
Smaller and Larger sizes can be manufactured upon request.



B10

Radshield's Lead Blockers are covered in our radiolucent range of non-tearable Vinyl or Nylon to ensure long life.

Available in:
.25mmPb, .35mmPb, .50mmPb or 1.0mmPb.

Standard Sizes:
24cm x 10cm
30cm x 12cm
43cm x 18cm

Note:
Smaller and Larger sizes can be manufactured upon request.

Lumbar Support Belt

Radshield's wide elasticised Lumbar Support Belt is designed for added comfort and support by reducing pressure off the back and shoulders and re-distributing the weight.

Standard Sizes:

XSmall – XLarge. Larger and Custom Made sizes available upon request.

Closure:

B1 = Adjustable Quick Release Buckle in the Front



B1

Mobile Height Adjustable Shields

Radshield's Mobile Height Adjustable Stands with Solid Lead Shields are ideal for patient and/or practitioner protection during diagnostic or image – amplification procedures, cath lab use and general protection. These Mobile Height Adjustable Shields move on a stable base with castors for stability and ease of rolling unit from one location to another. Height adjustable from 70cm to 162cm.

Available in:

.35mmPb, .50mmPb or 1.0mmPb.

Sizes Available:

Small: 61cm Wide x 76cm Long

Medium: 61cm Wide x 99cm Long

Large: 61cm Wide x 122cm Long

Note:

Larger and Custom Made sizes available upon request.



M4



M4 - Height Adjustable



M15



M15 - Height Adjustable

Curved Mobile Height Adjustable Trolley

Radshield's curved mobile height adjustable trolley, complete with .50mmPb lead shield. Height adjustment by means of springs and positioning pins. Complete with 4 castors (2 of which are lockable).

Available in:

.50mmPb

Sizes Available:

Minimum Height of Shield: 108cm

Maximum Height of Shield: 138cm

Opening Width: 55cm

Opening Depth: 53cm



M18



M18 - with lead shield

Ovarian Shields



O1

Radshield's Ovarian Shields are specifically for the protection of the genital area.

Available in:

.50mmPb or 1.0mmPb.

Set Comprises of:

- 2 x Small - 8cm x 14cm
- 1 x Medium - 15cm x 18cm
- 1 x Large - 18cm x 23cm

Note: Ovarian Shields can be purchased Individually or as a Set of 4.

Optional Extra's:

- Adjustable Hook-and-Loop Fastener Belt (**Code B3**)
- Wall-Mounted Hook-and-Loop Fastener Bar for out of the way storage (**Code VB**)



O2

Radshield's Ovarian Shields Set, complete with Adjustable Webbing Belt is specifically for the protection of the genital area.

Available in:

.50mmPb or 1.0mmPb.

Set Comprises of:

- 2 x Small - 8cm x 14cm
- 1 x Medium - 15cm x 18cm
- 1 x Large - 18cm x 23cm

Note: Ovarian Shields can be purchased Individually or as a Set of 4.

Optional Extra:

- Wall-Mounted Hook-and-Loop Fastener Bar for out of the way storage (**Code VB**).



O10

Radshield's Ovarian Shields complete with Ties are specifically for the protection of the genital area.

Available in:

.50mmPb or 1.0mmPb.

Set Comprises of:

- 2 x Small - 8cm x 14cm
- 1 x Medium - 15cm x 18cm
- 1 x Large - 18cm x 23cm

Note: Ovarian Shields can be purchased Individually or as a Set of 4.

Optional Extra:

- Wall-Mounted Hook-and-Loop Fastener Bar for out of the way storage (**Code VB**)



VB

Radshield's Hook-and-Loop Fastener Bar (for Wall-Mounting) for the storage of Gonad and Ovarian Shields.



B3

Radshield's fully Adjustable Webbing Belt is used with our Gonad (**Code G1**) and Ovarian (**Code O2**) Shields.

Closure:

B3 = Adjustable Quick Release Buckle

Paediatric Protection



C4-front



C4-back

Radshield's Paediatric Coat Type Wrap-Around apron's front design provides double protection for radiation sensitive sternal and gonadal areas.

Available in:

.25mmPb, .35mmPb or .50mmPb.

Standard Sizes:

Small: Age 1-3 (25cm Width x 66cm Length)

Medium: Age 4-6 (30cm Width x 71cm Length)

Large: Age 7-10 (36cm Width x 76cm Length)

Closure:

C4 = Adjustable Hook-and-Loop Fastener



S10-front



S10-back

Radshield's Paediatric Single-Sided Apron with Hook-and-Loop Fastener has always been a popular choice. That's because they're quick to put on, adjust and take off.

Available in:

.25mmPb, .35mmPb or .50mmPb.

Standard Sizes:

Small: Age 1-3 (30cm Width x 61cm Length)

Medium: Age 4-6 (36cm Width x 66cm Length)

Large: Age 7-10 (41cm Width x 71cm Length)

Closure:

S10 = Adjustable Hook-and-Loop Fastener at the Back
or

S5 = Adjustable Quick Release Buckle at the Back



M11

Radshield's Paediatric Mittens provide complete protection to the hands and are soft, supple and a pleasure to wear.

Available in:

25mmPb, .35mmPb or .50mmPb.

Standard Sizes:

One Size Fits All - Sold in Pairs



B7

Radshield's Paediatric Breast Shields protect baby's chest and lungs during spinal exams. Complete with Hook-and-Loop Fastener Neck Strap which allows for easy placement of the shields.

Available in:

.50mmPb.

Set Comprises of:

2 x Small 11cm x 14cm

2 x Medium 17cm x 20cm

2 x Large 19cm x 25cm



T5

Radshield's Paediatric Thyroid Collars provide protection to the thyroid and neck area.

Available in:

.50mmPb.

Standard Sizes:

Small: 8cm Width x 17cm Length

Medium: 9cm Width x 20cm Length

Large: 10cm Width x 23cm Length

Closure:

T5 = Adjustable Hook-and-Loop Fastener

Note: Thyroid Collars can be purchased Individually or as a Set of 3.



P3

Radshield's Paediatric Positioning Bags are covered in radiolucent material and double-stitched to prevent leakages. Each Positioning Bag includes a reinforced webbing handle.

Standard Sizes:

10cm x 13cm (.50kg)

10cm x 23cm (1kg)

13cm x 23cm (1.5kg)

18cm x 23cm (2kg)

15cm x 36cm (3kg)

Note:

Any size and/or shape can be manufactured upon request.

Paediatric Protection (Continued)



G8

Radshield's Paediatric Triangular Gonad Shields are specifically designed for the protection of the genital area.

Available in:
.50mmPb or 1.0mmPb.

Set Comprises of:
1 x Small 5cm x 5cm x 5cm
1 x Medium 7.5cm x 7.5cm x 7.5cm
1 x Large 10cm x 10cm x 10cm

Optional Extra's:

- Adjustable Hook-and-Loop Fastener Belt (**Code B3**)
- Wall-Mounted Hook-and-Loop Fastener Bar for out of the way storage of the Gonad Shields (**Code VBP**)



G10

Radshield's Paediatric Triangular Gonad Shields Set, complete with Adjustable Webbing Belt are specifically for the protection of the genital area.

Available in:
.50mmPb or 1.0mmPb.

Set Comprises of:
1 x Small - 5cm x 5cm x 5cm
1 x Medium - 7.5cm x 7.5cm x 7.5cm
1 x Large - 10cm x 10cm x 10cm

Note: Gonad Shields can be purchased Individually or as a Set of 3.

Optional Extra:

- Wall-Mounted Hook-and-Loop Fastener Bar for out of the way storage (**Code VBP**)



B3

Radshield's fully Adjustable Webbing Belt is used with our Gonad (**Code G1**) and Ovarian (**Code O2**) Shields.

Closure:
B3 = Adjustable Quick Release Buckle



VBP

Radshield's Paediatric Hook-and-Loop Fastener Bar (for Wall-Mounting) for the storage of Paediatric Gonad Shields (**Code G8**).



DC6-front



DC6-back

Radshield's Paediatric Dental Cape provides patients with optimum radiation protection.

Available in:
.25mmPb, .35mmPb or .50mmPb.

Standard Size:
Small: 50cm Wide x 50cm Long

Closure:
DC6 = Adjustable Hook-and-Loop Fastener

Note: Smaller and Custom Made sizes available upon request.

Patient Caliper

Radshield's lightweight aluminum Patient Caliper has an adjustable slide scale that reads in both inches and cm's.



C12

Patient Transfer Board

Transport your patients from bed to gurney to x-ray table safely and securely with **Radshield's** radiolucent Transfer Boards. The smooth surface slips easily under the patient and thanks to its radiolucent anti-static composition, you can leave it under the patient during the exam and scan. Suitable for use in CT, MRI or Radiology scans. Includes 4 (four) wrist straps.

Standard Size:

56cm Wide x 183cm Long x 0.50cm Thick



P4

Patient Slider, Complete with Padded Handles.

Transport your patients from bed to gurney to x-ray table safely and securely with **Radshield** Transfer Slider. Includes padded handles.

Standard Size:

72cm Wide x 105cm Long

Material:

Reinforced Oxford Nylon

Colour:

Black



PS1

Patient Wrap

Radshield's CT Patient Wrap protects the patient from unnecessary radiation exposure during CT studies. The Wrap goes under and around enveloping the patient. The convenient Hook-and-Loop Fastener and adjustable Webbing Straps provide easy and adjustable closure.

Available in:

.25mmPb, .35mmPb or .50mmPb

Standard Size:

155cm (Wide) x 76cm (Long)

Closure:

C2 = Adjustable Hook-and-Loop Fastener plus 2 x Webbing Straps with Quick Release Buckles

Note: Smaller and Larger Sizes can be manufactured upon request.



C2


Positioning Sponges



P2

Radshield's Covered Positioning Sponges are supplied in a Set of 10 assorted shapes and sizes for easy patient positioning. Covered in our radiolucent range of non-tearable Vinyl or Nylon to ensure long life.

Note: Positioning Sponges can be purchased Individually and any size and/or shape can be manufactured upon request.



P7

Radshield's Uncovered Positioning Sponges are supplied in a Set of 10 assorted shapes and sizes for easy patient positioning.

Note: Positioning Sponges can be purchased Individually and any size and/or shape can be manufactured upon request.

Positioning Bags



P3

Radshield's Positioning Bags are covered in our radiolucent range of non-tearable Vinyl, Nylon or PVC and are double-stitched to prevent any leakages. Each Positioning Bag includes a reinforced webbing handle for easy handling.

Standard Sizes:

- 10cm x 13cm (.50kg)
- 10cm x 23cm (1kg)
- 13cm x 23cm (1.5kg)
- 18cm x 23cm (2kg)
- 13cm x 30cm (2kg)
- 23cm x 23cm (2kg)
- 25cm x 25cm (3kg)
- 15cm x 36cm (3kg)
- 28cm x 28cm (4.5kg)
- 18cm x 41cm (4.5kg)
- 30cm x 30cm (4.5kg)
- 30cm x 36cm (5kg)
- 48cm x 30cm (6kg)

Note: Any size and/or shape can be manufactured upon request.



P6

Radshield's Positioning Bags are covered in our radiolucent range of non-tearable Vinyl, Nylon or PVC and are double-stitched to prevent any leakages. Each Positioning Bag includes a reinforced webbing handle for easy handling.

Standard Sizes:

- 10cm x 13cm (.50kg)
- 10cm x 23cm (1kg)
- 13cm x 23cm (1.5kg)
- 18cm x 23cm (2kg)
- 13cm x 30cm (2kg)
- 23cm x 23cm (2kg)
- 25cm x 25cm (3kg)
- 15cm x 36cm (3kg)
- 28cm x 28cm (4.5kg)
- 18cm x 41cm (4.5kg)
- 30cm x 30cm (4.5kg)
- 30cm x 36cm (5kg)
- 48cm x 30cm (6kg)

Note: Any size and/or shape can be manufactured upon request.

Protective Cap

Radshield's Protective Cap offers protection for your head against scatter radiation.

Available in:

.25mmPb, .35mmPb or .50mmPb.

Standard Sizes:

Large – XLarge. Larger and Custom Made sizes available upon request.

Note: For an accurate fit, please specify the circumference of your head just above your eyebrows.



S14-front



S14-back

Protective Shields/Blankets

Radshield's Protective Shields/Blankets are used for shielding patients and/or equipment.

Available in:

.25mmPb, .35mmPb, .50mmPb or 1.0mmPb.

Note: Any size and/or shape can be manufactured upon request.



P5

Scoliosis Stole

Radshield's Scoliosis Stole drapes over the shoulders with a webbing belt that will keep the Scoliosis Stole in place. The Scoliosis Stole provides breast and lung protection without compromising spinal image quality.

Available in:

.50mmPb.

Standard Sizes:

Large – XLarge. Larger and Custom Made sizes available upon request.



B2

Scrubs

Standard Sizes:

Small – XLarge. Larger and Custom Made sizes available upon request.

Note: Scrubs may be purchased as separate items. To order separately, kindly specify size and the piece you prefer. In addition, Top and/or Trouser sizes may be interchanged (i.e. small Top with large Trousers) at no additional charge. Please specify the size for each piece when ordering.



S9

Shin Guards

Radshield's Shin Guards give your tibia bone the ultimate protection. Secure them to your legs with adjustable Hook-and-Loop Fastener straps. Shin Guards are sold in pairs.

Available in:

.25mmPb, .35mmPb or .50mmPb.

Closure:

S8 = Adjustable Hook-and-Loop Fastener

Note: For an accurate fit, please specify calf circumference and preferred Shin Guard length.



S8

Storage

Hangers/Hooks



H2

Radshield's Wooden Hanger.



H3

Radshield's Metal Hanger.



H15

Radshield's stainless steel **Hanger Hooks** (set) can be fitted onto any of our steel hangers to correctly hang up Skirts, Half Aprons, etc.

Storage

Mobile Apron Trolleys



MT

Radshield's high quality **Mobile Apron Trolley**, complete with double-bars and storage tray, allows for storage of Aprons (bottom bar) and Thyroid Collars (top bar). The storage tray is used to store glasses, gloves etc. Designed with heavy-duty castors (2 of which are lockable), this unit can roll easily from one location to another and requires minimal space. Our **Steel (Code H3)** and **Wooden Hangers (Code H2)** work perfectly with this unit.

Standard Sizes:
190cm High x 95cm Wide x 55cm Deep



V10

Radshield's high quality Mobile Apron Trolley with heavy-duty castors (2 of which are lockable), is designed to roll easily from one location to another and requires minimal space. Our **Steel (Code H3)** and **Wooden (Code H2)** **Hangers** work perfectly with these units.

V10 = Holds 10 Aprons

Standard Sizes:
V10 - 189cm High x 57cm Wide x 50cm Deep



V11


Radshield's high quality **Mobile Apron Trolley** with heavy-duty castors (2 of which are lockable), is designed to roll easily from one location to another and requires minimal space. Our **Steel (Code H3)** and **Wooden (Code H2)** **Hangers** work perfectly with these units.

V11 = Holds 20 Aprons

Standard Sizes:
V11 - 189cm High x 115cm Wide x 50cm Deep

Storage

Wall-Mounted Racks



H4

Radshield's Single Apron Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.



H5

Radshield's Single Apron and 2 (Two) Glove Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.



H8

Radshield's 3 (Three) Apron and 2 (Two) Glove Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.



H13

Radshield's 2 (Two) Apron Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.




H23

Radshield's 2 (Two) Apron and 2 (Two) Glove Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.



H12

Radshield's 3 (Three) Apron Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.




H10

Radshield's 5 (Three) Apron Deluxe Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.




H14

Radshield's (Four) Apron Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.




H24

Radshield's 5 (Four) Apron Deluxe Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.



H7

Radshield's 5 (Five) Apron Economy Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.



H9

Radshield's 5 (Five) Apron Deluxe Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.



H16

Radshield's 6 (Six) Apron Economy Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.



H17

Radshield's 6 (Six) Apron Deluxe Wall-Mounted Rack, swings in either a left or right direction for out of the way storage.


Storage

Wall-Mounted Racks (Continued)




H28

Radshield's 1 (One) Apron Wall-Mounted Rack, complete with Hanging Adaptor, provides storage for Skirts and Half Aprons. Swings in either a left or right direction for out of the way storage.




H29

Radshield's 2 (Two) Apron Wall-Mounted Rack, complete with Hanging Adaptors, provides storage for Skirts and Half Aprons. Swings in either a left or right direction for out of the way storage.




H30

Radshield's 3 (Three) Apron Wall-Mounted Rack, complete with Hanging Adaptors, provides storage for Skirts and Half Aprons. Swings in either a left or right direction for out of the way storage.



H31

Radshield's 4 (Four) Apron Wall-Mounted Rack, complete with Hanging Adaptors, provides storage for Skirts and Half Aprons. Swings in either a left or right direction for out of the way storage.



H32

Radshield's 5 (Five) Apron Wall-Mounted Rack, complete with Hanging Adaptors, provides storage for Skirts and Half Aprons. Swings in either a left or right direction for out of the way storage.



H22

Radshield's Demi Half Apron Wall-Mounted Rack, holds up to 4 (Four) Demi Half Aprons.

Storage

Wall-Mounted Peg Racks



PR1
1-Apron Peg Rack
(Size: 200mm Long X 75mm Wide)



PR2
2-Apron Peg Rack
(Size: 300mm Long X 75mm Wide)



PR4
4-Apron Peg Rack
(Size: 508mm Long X 75mm Wide)



PR6
6-Apron Peg Rack
(Size: 767mm Long X 75mm Wide)



PR8
8-Apron Peg Rack
(Size: 1021mm Long X 75mm Wide)



PR10
10-Apron Peg Rack
(Size: 1275mm Long X 75mm Wide)



PR12
12-Apron Peg Rack
(Size: 1529mm Long X 75mm Wide)

Straps & Restraints



V5

Radshield's patient restraining Hook-and-Loop Fastener Straps are as strong as our Webbing Straps (**Code W5**) and soft on the skin. Easily applied, they are fully adjustable, re-useable and can be cut to any length.

Standard Sizes:

2.5cm Wide x 100cm Long
3.8cm Wide x 100cm Long
5cm Wide x 100cm Long
10cm Wide x 100cm Long

Note: Shorter and Longer lengths can be manufactured upon request.



W5

Radshield's patient restraining Webbing Straps are fully adjustable and re-useable.

Standard Sizes:

2.5cm Wide x 100cm Long
3.8cm Wide x 100cm Long
5cm Wide x 100cm Long

Note: Shorter and Longer lengths can be manufactured upon request.

Patient Restraint Belt



PRB

Radshield's patient restrain belt prevent retakes, unnecessary x-ray exposure and keep patients secure during x-ray examinations. Just wrap this belt around the patient to control his or her movements during the examination. The sturdy hook-and-loop fastener won't come undone and is strong enough to secure an unstable patient.

Standard Sizes:

13cm Wide x 195cm Long

Colour: Black

Note: Shorter and Longer lengths can be manufactured upon request.

Note: Shorter and Longer lengths can be manufactured upon request.

Table Drapes

Radshield's Table Drapes, which protect staff against scatter radiation to their lower extremities, can be manufactured to fit securely to any style of table. Select either a single-sided drape, or double the protection with double-sided drapes. When ordering, please specify the width and length you require.

Available in:

.25mmPb, .35mmPb, .50mmPb or 1.0mmPb.

Standard Sizes:

T6.1 91cm Wide x 61cm Long

T6.2 91cm Wide x 71cm Long

T6.3 91cm Wide x 91cm Long

T6.4 122cm Wide x 61cm Long

T6.5 122cm Wide x 71cm Long

T6.6 122cm Wide x 91cm Long

Note: Any size and/or shape can be manufactured upon request.



T6

Thyroid Collars

Radshield's Thyroid Collars offer protection to the thyroid, neck and sternum area.

Available in:

.50mmPb.

Standard Sizes:

Medium – 2XLarge



T1

T1 = Provides protection to the thyroid, neck and sternum area

Closure:

T1 = Hook-and-Loop Fastener at the Back

T1M = Magnet Closure



T3

T3 = Provides protection to the thyroid and neck area

Closure:

T3 = Hook-and-Loop Fastener at the Back

T3M = Magnet Closure



T13

T13 = Provides protection to the thyroid and neck area

Closure:

T13 = Hook-and-Loop Fastener at the Back

T13M = Magnet Closure



T7 = Provides protection to the thyroid, neck and sternum area

Closure:

T7 = Hook-and-Loop Fastener at the Back

T7M = Magnet Closure



T12

T12 = Provides protection to the thyroid, neck and sternum area

Closure:

T12 = Hook-and-Loop Fastener at the Back

T12M = Magnet Closure

Note: The T12 Thyroid Collar attaches to your Lead Apron.

Reusable Thyroid Collar Covers

Radshield's soft absorbent Thyroid Collar "Covers", protect Collars from perspiration and makeup.

Standard Sizes:

One Size Fits All

Note: Washable and Reusable.



T41

Fits over the **T1** Thyroid Collar



T43

Fits over the **T3** Thyroid Collar

Washable Aprons

Radshield's Reversible and Fully Washable Aprons

Standard Sizes:

Small – 60cm Wide
Medium – 70cm Wide
Large – 80cm Wide
XLarge – 90cm Wide

Note: Larger and Custom Made sizes available upon request.

Closure:

D12 = Ties at the Back



D12-front



D12-back

X-Ray Mattresses

Radshield's X-Ray Mattresses are covered in our radiolucent range of non-tearable Vinyl or Nylon to ensure long life. Our X-Ray Mattresses are manufactured with or without handles.

Standard Sizes:


185cm x 60cm x 3.5cm
220cm x 65cm x 3.5cm

Note:

Any size and/or shape can be manufactured upon request.



M10



Radshield

QUALITY X-RAY PROTECTIVE WEAR & ACCESSORIES

Contact Us

Telephone: 011 467 5577
011 467 5511

Email: radshield@mweb.co.za

Block C | Lonehill Office Park
Lone Close | Lonehill | Sandton



QUALITY X-RAY PROTECTIVE WEAR & ACCESSORIES